

# Low Elution 96- Well Magnet Plate

**All Permagen Labware Magnetic Separation Plates are designed for both performance and durability.**


- Solid aluminum construction for precise, trouble free use eliminating the concern for material inconsistencies or breakage associated with some plastic products on the market
- Precision engineered, giving you repeatable, reliable results
- Unique, fastener free construction, equals no screws or pins to come loose
- Side clearance for liquid handler grippers
- Hard-coat finished to maximize durability
- 24 Strong Neodymium magnets for fast separations of many magnetic bead chemistries  
\*beads are pulled to well sides
- SBS compliant footprint for automation
- Separations as low as 5  $\mu$ L
- Suitable for manual or robotic use
- Volume discounts available
- Most durable magnet plates on the market at competitive prices
- *Made in USA*

*For more information please visit our website or email us at [info@permagenlabware.com](mailto:info@permagenlabware.com)*

**Catalog no.: LE400**



**96- Well Low Elution Magnetic Separation Plate**



**Separations of 5  $\mu$ L - 50  $\mu$ L shown**



**Compatible with many 96- Well PCR plates**