



# 15 mL Centrifuge Tube Magnetic Separation Racks

**All Permagen Labware Centrifuge Tube Magnetic Separation Racks are designed for both performance and durability.**

- Solid aluminum alloy, one piece construction for trouble free use, eliminating the concern for material inconsistencies or breakage associated with some plastic products on the market
- Hard-coat finished to maximize durability
- Strong Neodymium magnets for fast separations of many magnetic bead chemistries
- Compatible with many 15 mL Centrifuge tubes
- Also available in configurations of:
  - (3) 50 mL Centrifuge, (3) 15 mL Centrifuge Combination Tube Rack
  - (1) 50 mL Centrifuge, (1) 15 mL Centrifuge Combination Tube Rack
  - (6) 50 mL Centrifuge Tube Rack
  - (2) 50 mL Centrifuge Tube Rack
- Volume discounts available
- Most durable separation racks on the market, at an affordable price
- *Made in USA*



**MSR6X15 – 6 X 15 mL Rack**



**MSR2X15 – 2 X 15 mL Rack**

*For more information please visit our website or email us at [info@permagenlabware.com](mailto:info@permagenlabware.com)*