



Microfuge Tube Magnetic Separation Racks

All Permagen Labware Microfuge Tube Magnetic Separation Racks are designed for both performance and durability.


- Solid, one piece aluminum construction eliminating the concern for material inconsistencies or breakage associated with some plastic products on the market for trouble free use
- Hard-coat finished to maximize durability
- Each tube has its own strong Neodymium magnet for fast separations of many magnetic bead chemistries
- Compatible with many 1.5 mL microfuge tubes and some 1.7 or 2 mL tubes
- Tube holes are angled directing 1.5 mL microfuge tubes to come in close proximity with the magnets
- Available in 6, 12, and 24 tube combinations
- Volume discounts available
- Most durable separation racks on the market, at an affordable price
- *Made in USA*




MSR06 – 6 X 1.5 mL Separation Rack



MSR12 – 12 x 1.5 mL Separation Rack



MSR24 – 24 X 1.5 mL Separation Rack



Magnetic Beads are pulled to inner tube wall

For more information please visit our website or email us at info@permagenlabware.com