

## *Permagen Magnetic Rack Line*

### *Separation Racks*

Permagen's separation racks are built to the same high standards as our magnet plates. They all feature a solid, one piece aluminum construction eliminating the concern of material failure associated with some racks currently on the market. All racks feature strong neodymium magnets. We offer many solutions for magnetic bead separations in Centrifuge Tubes, Microfuge Tubes, PCR Strips, and screw cap storage tubes.

Current offerings:

#### Centrifuge Tube Separation Racks

- 15 mL / 50 mL combination racks - Available in two configurations, one will hold up to (3) 15 mL and (3) 50 mL common Centrifuge tubes. The second option will hold (1) 15 mL and (1) 50 mL tube.
- 15 mL Racks - We offer two racks that will accommodate up to either (2 or 6) 15 mL Centrifuge Tubes.
- 50 mL Racks - Available in the same two configurations as above.

#### Microfuge Tube Separation Racks


- 1.5 mL Racks - Permagen offers three racks for magnetic bead separations in a 1.5 mL microfuge tube. They are available in a 6, 12, or 24 configuration.

#### PCR Strip Separation Rack

- 0.2 mL Strip - Will accommodate (2) 8 or 12 PCR strips with or without caps. Features angled tube holes assuring the tubes to be in close proximity to each magnet. Strong Neodymium magnets for fast separations.

#### Combination Rack

- 1.5 mL, PCR, 1.0 mL rack - This combination rack features 12 positions on one side which will allow separations in 1.5 mL microfuge tubes. The opposite side features 24 positions to accommodate either PCR strips, or 1.0 mL 2d style storage tubes.



For more information please visit our website  
or email us at [info@permagenlabware.com](mailto:info@permagenlabware.com)

Made in USA