

384- Post Magnet Plate


All Permagen Labware Magnetic Separation Plates are designed for both performance and durability.

- Solid aluminum construction for precise, trouble free use eliminating the concern for material inconsistencies or breakage associated with some plastic products on the market.
- Precision engineered, giving you repeatable, reliable results
- Fastener free construction, equals no screws or pins to come loose
- Side clearance for liquid handler grippers
- Hard-coat finished to maximize durability
- 96 Strong Neodymium post magnets for fast separations of many magnetic bead chemistries
*beads pulled to side of wells
- Compatible with many 384-Well Microplates
- Most durable magnet plates on the market at competitive prices
- SBS compliant footprint
- Manual or robotic use
- Maximum working volume of 39µl
- Volume discounts available
- *Made in USA*

For more information please visit our website or email us at info@permagenlabware.com



P384 – 384- Well Post magnet plate



Magnetic Beads are pulled to the sides of each well



Compatible with many 384- Well Microplates