


15 mL / 50 mL Centrifuge Tube Magnetic Separation Combo Racks


All Permagen Labware Centrifuge Tube Magnetic Separation Racks are designed for both performance and durability.

- Solid aluminum alloy, one piece construction for trouble free use, eliminating the concern for material inconsistencies or breakage associated with some plastic products on the market
- Hard-coat finished to maximize durability
- Strong Neodymium magnets for fast separations of many magnetic bead chemistries
- Compatible with many 15 mL and 50 mL Centrifuge tubes
- Also available in configurations of:
 - (6) 15 mL Centrifuge Tube Rack
 - (2) 15 mL Centrifuge Tube Rack
 - (6) 50 mL Centrifuge Tube Rack
 - (2) 50 mL Centrifuge Tube Rack
- Volume discounts available
- Most durable separation racks on the market, at an affordable price
- *Made in USA*

For more information please visit our website or email us at info@permagenlabware.com



MSR1550 – 3 X 15 mL/ 50 mL Rack



MSR50151 – 1 X 15 mL/ 50 mL Rack