

24- Tube Chiller Block

All Permagen Labware Products are designed for both performance and durability.

- Solid aluminum alloy construction for precise, trouble free use. Keep your samples chilled when in use or storage
- Precision engineered, giving you repeatable, reliable results
- Hard-coat finish to maximize durability
- Compatible with many 1.5 mL - 2.0 mL microfuge tubes and screw cap tubes
- SBS compliant footprint (127.7mm x 85.5mm)
- Manual or robotic use
- Fingertip groove around perimeter for secure grip
- Holes can be laser labeled alphanumerically as required
- Standard 18mm X 18mm tube offset for standard machine definitions
- Volume discounts available
- *Made in USA*

For more information please visit our website or email us at info@permagenlabware.com



Cb024 – 24- Tube Chiller Block



Compatible with many microfuge tubes and screw cap tubes

** shown with 2.0 mL microfuge tubes*