

M Series 96- Well Magnet Plate


All Permagen Labware Magnetic Separation Plates are designed for both performance and durability.

- Solid aluminum construction for precise, trouble free use eliminating the concern for material inconsistencies or breakage associated with some plastic products on the market.
- Precision engineered, giving you repeatable, reliable results
- Unique, fastener free construction, equals no screws or pins to come loose
- Side clearance for liquid handler grippers
- Hard-coat finished to maximize durability
- 96 Strong Neodymium ring magnets for fast separations of many magnetic bead chemistries
*beads are pulled into a small ring (approx. 4mm) near well bottoms
- SBS compliant footprint for automation
- Suitable for manual or robotic use with PCR plates and some microplates
- Available in two magnet strengths to suit your needs
-M480 – Up to 1 mL max or smaller volumes faster
-M380 – Up to 350 μ l
- Volume discounts available
- Most durable magnet plates on the market, at competitive prices
- *Made in USA*


For more information please visit our website or email us at info@permagenlabware.com



M series 96- Well Magnetic Separation Plate



Magnetic Beads are pulled into rings near the bottom of each well



Compatible with many 96- Well round or v-bottom Microplates and PCR plates